

MOB1A/MOB1B Rabbit Polyclonal Antibody



CAB18246

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

25kDa

Calculated MW:

Applications:

WB IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

55233/92597

Uniprot

Q9H8S9/Q7I9L4

Synonyms:

MOB1A/MOB1B


Immunogen:

Recombinant protein of human MOB1A/MOB1B.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - MOB1A/MOB1B Rabbit pAb (CAB18246)