

Vgll4 Rabbit Polyclonal Antibody



CAB18248

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

31KDa

Calculated MW:

31kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

232334

Uniprot

Q80V24

Synonyms:

VGL-4; BC048841; VGLL4


Immunogen:

Recombinant protein of human Vgll4.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - Vgll4 Rabbit pAb (CAB18248)