

CAB18248

Product Information

Product SKU:	CAB18248	Gene ID:	232334	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	31kDa	Conjugate:	-
Calculated MW:	31kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable transcription coactivator binding activity. Involved in several processes, including negative regulation of cardiac muscle cell proliferation; negative regulation of hippo signaling; and positive regulation of cellular protein catabolic process. Located in nucleus. Is expressed in lung; submandibular gland primordium; tooth; and vibrissa. Orthologous to human VGLL4 (vestigial like family member 4).

Recommended Dilution: WB,1:500 - 1:2000 IHC-P,1:50 - 1:200

Synonyms: VGL-4; Vgll4

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-287 of mouse Vgll4 (NP_808351.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.