

RMND1 Rabbit Polyclonal Antibody



CAB18254

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

51kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Protein Background

The protein encoded by this gene belongs to the evolutionary conserved sif2 family of proteins that share the DUF155 domain in common. This protein is thought to be localized in the mitochondria and involved in mitochondrial translation. Mutations in this gene are associated with combined oxidative phosphorylation deficiency-11. Alternatively spliced transcript variants have been found for this gene.

Immunogen information

Gene ID:

55005

Uniprot

Q9NWS8

Synonyms:

RMD1; C6orf96; COXPD11; bA351K16; bA351K16.3; RMND1

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

Recombinant protein of human RMND1.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
