

CAB18260

Product Information

Product SKU:	CAB18260	Gene ID:	54468	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	100-102kDa	Conjugate:	-
Calculated MW:	99kDa	Isotype:	IgG

Immunogen Information

Background: Involved in cellular response to amino acid starvation; positive regulation of TOR signaling; and protein-containing complex localization. Located in several cellular components, including cytosol; lysosomal membrane; and nucleoplasm. Part of GATOR2 complex.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: MIO; Sea4; Yulink; MIOS

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 556-875 of human MIOS (NP_061878.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.