

CAB18261

Product Information

Product SKU:	CAB18261	Gene ID:	-	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat,Other (Wide Range Predicted)

Additional Information

Observed MW:	22-70kDa	Conjugate:	-
Calculated MW:	-	Isotype:	IgG

Immunogen Information

**Background:** The methylation of arginine residues in proteins is a post-translational modification that contributes to a wide range of biological processes. Many cytokines involved in T cell development and activation utilize the common cytokine receptor  $\gamma$ -chain ( $\gamma c$ ) and the kinase JAK3 for signal transduction.

Recommended Dilution:	WB,1:100 - 1:500
Synonyms:	-
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence containing succinylated R.
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.