

PLXNA1 Rabbit Polyclonal Antibody



CAB18264

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

211kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

5361

Uniprot

Q9UIW2

Synonyms:

NOV; NOV; PLXN1; PLEXIN-A1; PLXNA1

Immunogen:

A synthetic peptide of human PLXNA1.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
