

# GSDMD Rabbit Polyclonal Antibody



CAB18281

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

53KDa/23KDa

### Calculated MW:

53kDa

### Applications:

WB IHC IF IP

### Reactivity:

Human, Mouse, Rat

## Protein Background

Gasdermin D is a member of the gasdermin family. Members of this family appear to play a role in regulation of epithelial proliferation. Gasdermin D has been suggested to act as a tumor suppressor. Alternatively spliced transcript variants have been described.

## Immunogen information

### Gene ID:

79792

### Uniprot

P57764

### Synonyms:

DF5L; DFNA5L; FKSG10; GSDMDC1; GSDMD

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50  
- 1:200 IF 1:50 - 1:200 IP  
1:50 - 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

### Immunogen:


Recombinant protein of human GSDMD.

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot - GSDMD Rabbit pAb (CAB18281)