

LGALS7 Rabbit Polyclonal Antibody



CAB18284

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

15kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Protein Background

The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. Differential and in situ hybridization studies indicate that this lectin is specifically expressed in keratinocytes and found mainly in stratified squamous epithelium. A duplicate copy of this gene (GenelD:653499) is found adjacent to, but on the opposite strand on chromosome 19.

Immunogen information

Gene ID:

3963

Uniprot

P47929

Synonyms:

GAL7; LGALS7A; LGALS7

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

Recombinant protein of human LGALS7.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
