

Product Information

Product SKU:	CAB18321	Gene ID:	8799	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	28kDa	Conjugate:	-
Calculated MW:	28kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene facilitates peroxisomal proliferation and interacts with PEX19. The encoded protein is found in the peroxisomal membrane. Several transcript variants, some protein-coding and some not protein-coding, have been found for this gene.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	PEX14B; PEX11beta; PEX11-BETA; PEX11B
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 100-180 of human PEX11B (NP_003837.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.