

# TPX2 Rabbit Polyclonal Antibody



CAB18327

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

105kDa

### Calculated MW:

85kDa

### Applications:

WB IF

### Reactivity:

Human, Rat

## Protein Background

### Immunogen information

#### Gene ID:

22974

#### Uniprot

Q9ULW0

#### Synonyms:

DIL2; p100; DIL-2; HCTP4; FLS353; HCA519; REPP86; C20orf1; C20orf2; GD:C20orf1; TPX2

#### Immunogen:

Recombinant protein of human TPX2.

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:


Affinity purification

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot - TPX2 Rabbit pAb (CAB18327)