

WLS Rabbit Polyclonal Antibody



CAB18328

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

70kDa

Calculated MW:

62kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

79971

Uniprot

Q5T9L3

Synonyms:

EVI; MRP; GPR177; mig-14; C1orf139; WLS


Immunogen:

Recombinant protein of human WLS.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - WLS Polyclonal Antibody (CAB18328)