

MLLT1 Rabbit Polyclonal Antibody



CAB18333

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

80kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

4298

Uniprot

Q03111

Synonyms:

ENL; LTG19; YEATS1; MLLT1


Immunogen:

Recombinant protein of human MLLT1.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - MLLT1 Polyclonal Antibody (CAB18333)