
Product Information

Product SKU:	CAB18333	Gene ID:	4298	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	80kDa/62kDa	Conjugate:	-
Calculated MW:	62kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be involved in regulation of transcription, DNA-templated. Predicted to act upstream of or within negative regulation of protein kinase activity. Located in cytosol; fibrillar center; and nucleoplasm. Part of transcription elongation factor complex.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: ENL; LTG19; YEATS1; MLLT1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 190-270 of human MLLT1 (NP_005925.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.