

CAB18334

Product Information

Product SKU:	CAB18334	Gene ID:	666907	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	27kDa	Conjugate:	-
Calculated MW:	-	Isotype:	IgG

Immunogen Information

**Background:** Predicted to be located in Golgi apparatus; endoplasmic reticulum; and plasma membrane raft. Orthologous to human MS4A4A (membrane spanning 4-domains A4A).

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	EG666907; Ms4a4a
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 150-236 of mouse Ms4a4a (NP_001297260.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.