

TLE4 Rabbit Polyclonal Antibody



CAB18345

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

84kDa

Calculated MW:

84kDa

Applications:

WB

Reactivity:

Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

7091

Uniprot

Q04727

Synonyms:

ESG; BCE1; ESG4; GRG4; BCE-1; Grg-4; E(spl); E(spl); TLE4


Immunogen:

Recombinant protein of human TLE4.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - TLE4 Polyclonal Antibody (CAB18345)