

RECQL5 Rabbit Polyclonal Antibody



CAB18346

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

9400

Uniprot

O94762

Synonyms:

ESG; BCE1; ESG4; GRG4; BCE-1; Grg-4; E(spl); E(spl); RECQL5; RECQ5

Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 400-500 of human RECQL5 (NP_004250.4).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
