

XRN2 Rabbit Polyclonal Antibody



CAB18350

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

108kDa

Calculated MW:

108kDa

Applications:

WB

Reactivity:

Human

Protein Background

This gene encodes a 5'-3' exonuclease that promotes transcription termination at cotranscriptional cleavage sites. Alternative splicing results in multiple transcript variants encoding different isoforms.

Immunogen information

Gene ID:

22803

Uniprot

Q9H0D6

Synonyms:

XRN2

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification


Immunogen:

Recombinant protein of human XRN2.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - XRN2 Rabbit pAb (CAB18350)