

CELSR2 Rabbit Polyclonal Antibody



CAB18359

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

WB

Reactivity:

Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

1952

Uniprot

Q9HCU4

Synonyms:

ADGRC2; CDHF10; EGFL2; Flamingo1; MEGF3; CELSR2

Immunogen:

Recombinant protein of human CELSR2.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
