

HOXD13 Rabbit Polyclonal Antibody



CAB18365

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

3239

Uniprot

P35453

Synonyms:

BDE; BDSL; HOX4I; SPD; SPD1; HOXD13

Immunogen:

Recombinant protein of human HOXD13.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
