

# Polyclonal AntibodyPC3 Rabbit Polyclonal Antibody

## CAB18370



### Product Information

**Size:**

20uL, 50uL, 100uL, 200uL

**Observed MW:**

Refer to figures

**Calculated MW:**

70kDa

**Applications:**

WB

**Reactivity:**

Human, Mouse

### Antibody Information

**Recommended dilutions:**

WB 1:500 - 1:2000

**Source:**

Rabbit

**Isotype:**

IgG

**Purification:**

Affinity purification

### Protein Background

Messenger RNA stability and translation initiation are extensively under the control of poly(A)-binding proteins (PABP). See PABPC1 (MIM 604679) for background information.

### Immunogen information

**Gene ID:**

5042

**Uniprot**

Q9H361

**Synonyms:**

PABP3; PABPL3; tPABP; PABPC3

**Immunogen:**

Recombinant protein of human PABPC3.

**Storage:**

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---