

CAB18373

Product Information

Product SKU:	CAB18373	Gene ID:	5586	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	140kDa	Conjugate:	-
Calculated MW:	112kDa	Isotype:	IgG

Immunogen Information

Background: Enables RNA polymerase binding activity; histone deacetylase binding activity; and protein serine/threonine kinase activity. Involved in several processes, including apical junction assembly; positive regulation of cell cycle; and positive regulation of viral genome replication. Located in several cellular components, including cleavage furrow; cytoskeleton; and midbody.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: PAK2; PRK2; STK7; Pak-2; PRKCL2; PRO2042; PRK2/PKN2

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 600-700 of human PRK2/PRK2/PKN2 (NP_006247.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.