

PKN2 Rabbit Polyclonal Antibody



CAB18373

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

140kDa

Calculated MW:

75kDa/94kDa/106kDa/110kDa/112kDa

Applications:

WB

Reactivity:

Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

5586

Uniprot

Q16513

Synonyms:

PAK2; PRK2; PRKCL2; PRO2042; Pak-2; STK7; PKN2


Immunogen:

Recombinant protein of human PKN2.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - PKN2 Rabbit pAb (CAB18373)