

SNTB2 Rabbit Polyclonal Antibody



CAB18384

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

6645

Uniprot

Q13425

Synonyms:

D16S2531E; EST25263; SNT2B2; SNT3; SNTL; SNTB2

Immunogen:

Recombinant protein of human SNTB2.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
