

CAB18400

Product Information

Product SKU:	CAB18400	Gene ID:	8843	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	44kDa	Conjugate:	-
Calculated MW:	44kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable GTP binding activity and purinergic nucleotide receptor activity. Predicted to be involved in G protein-coupled receptor signaling pathway. Located in cell junction.

Recommended Dilution:	WB,1:500 - 1:2000 IHC-P,1:50 - 1:200
Synonyms:	HCA3; HM74; PUMAG; Puma-g; GPR109B; HCAR3
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 300-392 of human HCAR3 (NP_006009.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.