

HCAR3 Rabbit Polyclonal Antibody



CAB18400

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

44kDa

Calculated MW:

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

8843

Uniprot

P49019

Synonyms:

GPR109B; HCA3; HM74; PUMAG; Puma-g; HCAR3

Immunogen:

Recombinant protein of human HCAR3.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot - HCAR3 Rabbit pAb (CAB18400)

