

# TRIM66 Rabbit Polyclonal Antibody



CAB18411

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

135kDa

### Calculated MW:

### Applications:

WB

### Reactivity:

Human

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

9866

#### Uniprot

O15016

#### Synonyms:

C11orf29; TIF1D; TIF1DELTA; TRIM66

#### Immunogen:

Recombinant protein of human TRIM66.

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---

Western blot - TRIM66 Rabbit pAb (CAB18411)

