

Product Information

Product SKU:	CAB18411	Gene ID:	9866	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	135kDa	Conjugate:	-
Calculated MW:	149kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable chromatin binding activity and identical protein binding activity. Predicted to act upstream of or within negative regulation of transcription, DNA-templated. Located in aggresome and nucleoplasm.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	TIF1D; C11orf29; TIF1DELTA; TRIM66
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 750-970 of human TRIM66 (NP_055633.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.