

SH2D3C Rabbit Polyclonal Antibody



CAB18416

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

10044

Uniprot

Q8N5H7

Synonyms:

CHAT; NSP3; PRO34088; SHEP1; SH2D3C

Immunogen:

Recombinant protein of human SH2D3C.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
