

# B3GNT2 Rabbit Polyclonal Antibody



CAB18418

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

### Applications:

WB

### Reactivity:

Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

10678

#### Uniprot

Q9NY97

#### Synonyms:

3-Gn-T1; 3-Gn-T2; B3GN-T2; B3GNT; B3GNT-2; B3GNT1;  
BETA3GNT; BGNT2; BGnT-2; beta-1; beta3Gn-T1; beta3Gn-T2;  
B3GNT2

#### Immunogen:

Recombinant protein of human B3GNT2.

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---