

C1QL1 Rabbit Polyclonal Antibody



CAB18419

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

26kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

10882

Uniprot

O75973

Synonyms:

C1QRF; C1QTNF14; CRF; C1QL1


Immunogen:

A synthetic peptide of human C1QL1.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - C1QL1 Polyclonal Antibody (CAB18419)