

Product Information

Product SKU:	CAB18419	Gene ID:	10882	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	26kDa	Conjugate:	-
Calculated MW:	26kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable signaling receptor binding activity. Predicted to act upstream of or within maintenance of synapse structure; motor learning; and neuron remodeling. Predicted to be located in several cellular components, including climbing fiber; presynapse; and synaptic cleft. Predicted to be part of collagen trimer.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	CRF; C1QRF; CTRP14; C1QTNF14; C1QL1
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 38-88 of human C1QL1 (NP_006679.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.