

SHANK2 Rabbit Polyclonal Antibody



CAB18423

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

159kDa

Calculated MW:

Applications:

WB IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

22941

Uniprot

Q9UPX8

Synonyms:

AUTS17; CORTBP1; CTTNBP1; ProSAP1; SHANK; SPANK-3; SHANK2


Immunogen:

Recombinant protein of human SHANK2.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - SHANK2 Rabbit pAb (CAB18423)