
Product Information

Product SKU:	CAB18427	Gene ID:	23165	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	250kDa	Conjugate:	-
Calculated MW:	228kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a nucleoporin, which is a subunit of the nuclear pore complex that functions in active transport of proteins, RNAs and ribonucleoprotein particles between the nucleus and cytoplasm. Mutations in this gene are associated with steroid-resistant nephrotic syndrome.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: NPHS13; C7orf14; NUP205

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 10-120 of human NUP205 (NP_055950.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.