

ACSBG1 Rabbit Polyclonal Antibody



CAB18431

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

23205

Uniprot

Q96GR2

Synonyms:

BG; BG1; BGM; GR-LACS; LPD; ACSBG1

Immunogen:

Recombinant protein of human ACSBG1.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
