

---

## Product Information

<b>Product SKU:</b>	CAB18432	<b>Gene ID:</b>	23221	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

## Additional Information

<b>Observed MW:</b>	83kDa	<b>Conjugate:</b>	-
<b>Calculated MW:</b>	83kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** The protein encoded by this gene is a small Rho GTPase and a candidate tumor suppressor. The encoded protein interacts with the cullin-3 protein, a ubiquitin E3 ligase necessary for mitotic cell division. This protein inhibits the growth and spread of some types of breast cancer. Three transcript variants encoding different isoforms have been found for this gene.

**Recommended Dilution:** WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200

**Synonyms:** p83; DBC2; DEE64; EIEE64; RHOBTB2

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 320-480 of human RHOBTB2 (NP\_001153508.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.