

FAM32A Rabbit Polyclonal Antibody



CAB18444

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

WB

Reactivity:

Human

Protein Background

Immunogen information

Gene ID:

26017

Uniprot

Q9Y421

Synonyms:

OTAG-12; OTAG12; FAM32A

Immunogen:

Recombinant protein of human FAM32A.

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
