

SIPA1L1 Rabbit Polyclonal Antibody
CAB18445



Product Information

Product SKU:	CAB18445	Gene ID:	26037	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	180kDa	Conjugate:	-
Calculated MW:	200kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable GTPase activator activity; actin filament binding activity; and protein kinase binding activity. Predicted to be involved in several processes, including actin cytoskeleton organization; activation of GTPase activity; and regulation of postsynapse organization. Located in actin cytoskeleton and plasma membrane.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	E6TP1; SPAR1; SIPA1L1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1480-1700 of human SIPA1L1 (NP_056371.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.