

Product Information

Product SKU:	CAB18451	Gene ID:	29803	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	64kDa	Conjugate:	-
Calculated MW:	64kDa	Isotype:	IgG

Immunogen Information

Background:	Enables RNA binding activity. Predicted to be involved in regulation of transcription by RNA polymerase II. Predicted to act upstream of or within positive regulation of glucose import and regulation of fatty acid transport. Located in nucleoplasm.
Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	AP4; RIP60; ZNF464; Zfp464; REPIN1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 376-567 of human REPIN1 (NP_037532.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.