

Product Information

Product SKU:	CAB18454	Gene ID:	51109	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	35kDa	Conjugate:	-
Calculated MW:	35kDa	Isotype:	IgG

Immunogen Information

Background: Enables NADP-retinol dehydrogenase activity. Involved in cellular detoxification of aldehyde and retinoid metabolic process. Acts upstream of or within retinal metabolic process. Located in endoplasmic reticulum membrane.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	MDT1; CGI82; PSDR1; RALR1; SCALD; ARSDR1; HCBP12; RDJCSS; SDR7C1; RDH11
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 70-145 of human RDH11 (NP_057110.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.