

# UTP11L Rabbit Polyclonal Antibody



CAB18455

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

30kDa

### Calculated MW:

### Applications:

WB

### Reactivity:

Human, Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

51118

#### Uniprot

Q9Y3A2

#### Synonyms:

CGI-94; CGI94; UTP11L

#### Immunogen:

Recombinant protein of human UTP11L.

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---

Western blot - UTP11L Rabbit pAb (CAB18455)

