

Product Information

Product SKU:	CAB18460	Gene ID:	51490	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	42kDa	Conjugate:	-
Calculated MW:	42kDa	Isotype:	IgG

Immunogen Information

Background: Enables miRNA binding activity. Involved in maintenance of centrosome location and production of miRNAs involved in gene silencing by miRNA. Located in kinetochore; mitotic spindle; and spindle pole centrosome.

Recommended Dilution:	WB,1:500 - 1:2000 IHC-P,1:50 - 1:200
Synonyms:	CENP32; CENP-32; HSPC109; C9orf114
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 163-340 of human C9orf114 (NP_057474.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.