

C9orf114 Rabbit Polyclonal Antibody



CAB18460

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

42kDa

Calculated MW:

Applications:

WB IHC

Reactivity:

Human, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

51490

Uniprot

Q5T280

Synonyms:

C9orf114; CENP-32; HSPC109


Immunogen:

Recombinant protein of human C9orf114.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - C9orf114 Rabbit pAb (CAB18460)