

UQCC1 Rabbit Polyclonal Antibody



CAB18469

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

WB

Reactivity:

Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

55245

Uniprot

Q9NVA1

Synonyms:

BFZB; C20orf44; CBP3; UQCC; UQCC1

Immunogen:

Recombinant protein of human UQCC1.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
