

# UNC45A Rabbit Polyclonal Antibody



CAB18474

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

### Applications:

WB

### Reactivity:

Human

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

55898

#### Uniprot

Q9H3U1

#### Synonyms:

GC-UNC45; GCUNC-45; GCUNC45; IRO039700; SMAP-1; SMAP1; UNC-45A; UNC45A

#### Immunogen:

Recombinant protein of human UNC45A.

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---