

Product Information

Product SKU:	CAB18475	Gene ID:	55900	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	55kDa	Conjugate:	-
Calculated MW:	55kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the zinc-finger protein family. The encoded protein contains seven C2H2-type zinc fingers and a KRAB domain, but its function has yet to be determined. Alternatively spliced transcript variants have been described.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	HSD16; MST154; ZNF327; MSTP154; ZNF135L; ZNF140L; ZNF302
Purification Method:	Affinity purification
Immunogen:	Recombinant Protein corresponding to a sequence within amino acids 90-200 of human ZNF302 (NP_060913.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.