

# UBL5 Rabbit Polyclonal Antibody



CAB18481

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

### Applications:

IF

### Reactivity:

Human, Mouse, Rat

## Antibody Information

### Recommended dilutions:

IF 1:50 - 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

59286

#### Uniprot

Q9BZL1

#### Synonyms:

HUB1; UBL5

#### Immunogen:

Recombinant protein of human UBL5.

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---