

LMBR1 Rabbit Polyclonal Antibody



CAB18484

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

64327

Uniprot

Q8WVP7

Synonyms:

ACHP; C7orf2; DIF14; LSS; PPD2; THYP; TPT; ZRS; LMBR1

Immunogen:

Recombinant protein of human LMBR1.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
