

# Rptor Rabbit Polyclonal Antibody



CAB18487

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

## Protein Background

### Immunogen information

#### Gene ID:

74370

#### Uniprot

### Applications:

WB

### Reactivity:

Human

### Synonyms:

4932417H02Rik; Rap; Raptor; mKIAA1303; Rptor; RPTOR; KOG1; Mip1

### Immunogen:

Recombinant protein of human Rptor.

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---