

CHCHD7 Rabbit Polyclonal Antibody



CAB18489

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

10kDa

Calculated MW:

6kDa/7kDa/10kDa/11kDa/13kDa

Applications:

WB IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

CHCHD7 (coiled-coil-helix-coiled-coil-helix domain containing 7) is a protein-coding gene.

Immunogen information

Gene ID:

79145

Uniprot

Q9BUK0

Synonyms:

COX23; CHCHD7

Immunogen:

Recombinant protein of human CHCHD7.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot - CHCHD7 Rabbit pAb (CAB18489)

